

## DATA SHEET

### Standard Vicat Apparatus



63-L0028

#### Descripción General

The Vicat frame consists essentially of a metal stand with a sliding rod. An adjustable indicator moves over a graduated scale. The needle or plunger is attached to the bottom end of the rod to make up the test weight of 300 g.

The frame 63-L0028/1 is supplied without accessories, which have to be ordered separately depending on the requirement. It is also proposed with basic EN accessories (63-L028) or ASTM/AASHTO accessories (63-L0028/A).

- Weight approx.: 4 kg

Note: For the automatic version of the Vicat apparatus see Vicamatic

#### Normativas

- EN 196-3
- EN 480-2
- ASTM C191
- AASHTO T131

## DATA SHEET



Accessories

### Productos

#### 63-L0028

Vicat test set EN method

#### 63-L0028/A

Vicat test set ASTM method

#### 63-L0028/1

Vicat apparatus

### Accesorios y Consumibles

#### 63-L0028/E18

Accessory for Vicat apparatus series L0028 comprehending conical penetration probe 8 mm dia x50 mm long with 100g calibrated weight for gypsum testing to EN 13279

#### 63-L0028/21

Initial/final needle 1 mm dia set of 6

#### 63-L0028/31

Initial needle 1,13 mm dia set of 6

#### 63-L0028/41

Stainless steel final needle with special foot (single piece)

#### 63-L0028/5

Consistency plunger 10 mm dia

#### 63-L0028/6

Supporting glass plate

#### 63-L0028/7

Glass thermometer -10 + 50c

#### 63-L0028/8

Additional weight 700 g

#### 63-L0027/E22

## DATA SHEET

---

Plastic mould to EN

**63-L0027/E23**

Plastic mould to ASTM/AASHTO